






Generation **3**



The Orthocare Range of Orthopaedic



















About Us



Accordingly, ORTHO CARE, a USFDA 510K Cleared & CE certify company with manufacturing practices and facilities that conform to the guidelines of ISO 13485, is in the business of manufacturing of Orthopedic Implants. Also The company produces an array of quality implants which are affordable and are trusted by healthcare professionals and patients across geographies. However, In a very short span of time, as can be seen being a leader in the domestic market, though, the company has marked global footprints in 60 countries across the world. As a matter of fact, The Company is focusing on increasing the momentum in its respective category of business in the international market at an aggressive pace.

If you are looking for power tools that offer reliable performance in the most demanding applications, that come with a variety of options to suit your needs and are available at reasonable prices, ORTHO CARE is the answer.



Generation 3

1

2

3

4

5

6

7

8

9-10

11-12

13-15

Ergonomical Shape

Index

1. Cannulated Di	rill Machine

2. Oscillating Saw Machine

3. Sternum Saw/Reciprocating Saw

4. Acetabular Reaming Drill

5. Multifunction Ortho Dril & Saw

6. Micro Cannulated Drill

7. Micro Oscillating Saw

8. Micro Multifunction Drill & Saw

9. Saw Blades 10. K-Wire Attachment

11. Spares

Long Battery life

Faster Speed

High Power

Excellent Grip



Cannulated Drill Machine









With use of G3.217 attachment (option)



With use of G3.219 attachment (option)

It is used for various types of drilling operations. Also suitable for K-wire surgery *& Reaming Operations*

Standard Configuration:

⋆ One Hand Piece

★ 1 pc. Charger

★ 2 pcs. Battery

★ 2 pcs. Sterilizing Rings

* 1 pc. chuck & key

Technical Specifications:

	poomoations.		
Speed	1000rpm		
Noise	≤65db		
Sterilize	Fully Autoclavable (without battery) at 121°C/135°C		
Torque	≥3.5 NM		
Han d	Working Voltage	14.4V	
Hand Piece	Output Power	≥90W	
i iece	Unit Rise in Temp.	≤20°C	
Charger	Input Voltage	AC 110V-240V 50HZ-60HZ	
Charge	Output Power	7.2V - 14.4V	
	Voltage	14.4V	
Battery	Charging Time 2Hours or more		
	Running Time	45-50 minutes on full charge	
Packing	Packed in Elega	int Carrying Case	

^{*} With use of optional attachments

Battery Operated







Standard Configuration:

★ One Hand Piece

★1 pc. Charger

★ 2 pcs. Saw Blades

★2 pcs. Sterilizing Rings

★ 2 pcs. Battery

Speed	0-20000 osc/r	n	
Noise	≤70db		
Sterilize	Fully Autoclavable (without battery) at 121°C/135°C		
Hand	Working Voltage	14.4V	
Hand Piece	Output Power	≥90W	
1 lece	Unit Rise in Temp.	≤20°C	
Charger	Input Voltage	AC 110V-240V 50HZ-60HZ	
Charger	Output Power	14.4V	
	Voltage	14.4V	
Battery	Charging Time	2Hours or more	
	Running Time 45-50 minutes on full charge		
Packing	Packed in Elegant Carrying Case		





Sternum Saw/Reciprocating Saw





Standard Configuration:

- * One Hand Piece
- * Sternum Top
- * 2 pcs. Sternum Saw Blades
- ★ 2 pcs. Reciprocating Saw Blades
- * 1 pc. Charger
- * 2 pcs. Sterilizing Rings
- * 2 pcs. Battery

Speed	0-15000 osc/r	n	
Noise	≤75db		
Sterilize	Fully Autoclavable (without battery) at 121°C/135°C		
Hand	Working Voltage	14.4V	
Piece	Output Power	≥90W	
1 1000	Unit Rise in Temp.	≤20°C	
Charger	Input Voltage	AC 110V-240V 50HZ-60HZ	
Citarger	Output Power	7.2V - 14.4V	
	Voltage	14.4V	
Battery	Charging Time	2Hours or more	
	Running Time 45-50 minutes on full charge		
Packing	Packed in Elegant Carrying Case		





Acetabular Reaming Drill





Acetabular reaming drill can be used in total hip artlopplasty to prepare a patient's acetabulum for implantation of a hip prosthesis

Standard Configuration:

★ One Hand Piece

★ 1 pc. Charger

★2 pcs. Battery

★ 2 pcs. Sterilizing Rings

*1 pc. chuck & key

Technical Specifications:

	poomoations.		
Speed	300rpm		
Noise	≤65db		
Sterilize	Fully Autoclav	able (without battery) at 121°C/135°C	
Torque	≥6.5 NM		
Hand Piece	Working Voltage Output Power Unit Rise in Temp.	14.4V ≥90W <20°C	
	Input Voltage	AC 110V-240V 50HZ-60HZ	
Charger	Output Power	7.2V - 14.4V	
Battery	Voltage Charging Time Running Time	14.4V 2Hours or more 45-50 minutes on full charge	
Packing	Packing Packed in Elegant Carrying Case		

^{*} Reamers not included.



Ergonomical Shape



Multifunction Ortho Drill & Saw





Standard Configuration:

- ★ Hand Piece Motor, Battery Operated
- ★ Drill Attachment (Jacob Chuck Tpe)
- ★ Cannualted Drill Attachment
- ♦ Oscillating Saw attachment
- ★ Reamer Attachment
- * Convert Adapter for Quick Coupling Drill Bits
- * Convert Adapter for AO design Reamers
- ★1 pc. Charger
- ★ 2 pcs. Sterilizing Rings
- ★ 2 pcs. Oscillating Saw Blades
- ★ 2 pcs. Battery

Technical Specifications:

Speed Drill	800rpm		
Speed Saw	20,000ocs/m		
Noise	≤65db		
Sterilize	Fully Autoclavable (without battery) at 121°C/135°C		
Torque	≥3.5 NM		
Hand Piece	Working Voltage Output Power Unit Rise in Temp.	14.4V ≥90W <20°C	
Charger	Input Voltage Output Power	AC 110V-240V 50HZ-60HZ 7.2V - 14.4V	
Battery	Voltage Charging Time Running Time	14.4V 2Hours or more 45-50 minutes on full charge	
Packing	Packed in Elega	nt Carrying Case	



Battery Operated

Ergonomical Shape



Micro Cannulated Drill





Standard Configuration:

- ★ One hand piece
- ★ Two battery7.2V
- ★ One charger

- ★ One drill chuck and 1 key
- ★ Two Sterilizing Rings
- ★ One aluminum box

Technical Specifications:

Speed	1000rpm		
Noise	≤63db		
Sterilize	Fully Autoclav	able (without battery) at 121°C/135°C	
Torque	≥3.3 NM		
Hand Piece	Working Voltage Output Power Unit Rise in Temp.	7.2V ≥50W ≤20°C	
Charger	Input Voltage Output Power	AC 110V-220V 50HZ-60HZ 7.2V - 14.4V	
Battery	Voltage Charging Time Running Time	7.2V 2Hours or more 45-50 minutes on full charge	
Packing	Packed in Elegant Carrying Case		

Ergonomical Shape



Micro Oscillating Saw





Standard Configuration:

- ★ One hand piece
- ★ Two battery 7.2V
- ★ One charger

- ★ Two saw blades
- ★ Two Sterilizing Rings
- ★ One aluminum box

Speed	0-18000 osc/r	n	
Noise	≤65db		
Sterilize	Fully Autoclav	able (without battery) at 121°C/135°C	
Hand Piece	Working Voltage 7.2V Output Power ≥50W Unit Rise in Temp. ≤20°C		
Charger	Input Voltage Output Power	AC 110V-220V 50HZ-60HZ 7.2V - 14.4V	
Battery	Voltage 7.2V Charging Time 2Hours or more Running Time 45-50 minutes on full charge		
Packing	Packed in Elegant Carrying Case		





Micro Multifunction Drill & Saw G3 - 503

Standard Configuration:

- ★ Cannulated Hand Piece Motor, Battery Operated
- ★ Drill attachment (with Jacob Chuck)
- * Convert Adapter for Quick Coupling Drill bits
- ★ K-Wire attachment for wires 0.8mm- 2.5mm
- ★ Oscillating Saw attachment

- ★ Oscillating Saw attachment
- ★ 5 oscillating Saw Blades
- ★2 pcs. Battery
- ★1 pc. Charger
- * 2 pcs Sterilizing Ring

Speed Drill	800rpm		
Speed Saw	18,000ocs/m		
Noise	≤65db		
Sterilize	Fully Autoclav	able (without battery) at 121°C/135°C	
Torque	≥3.3 NM		
111	Working Voltage	7.2V	
Hand Piece	Output Power	≥50W	
riece	Unit Rise in Temp.	≤20°C	
Charger	Input Voltage AC 110V-220V 50HZ-60HZ		
Charger	Output Power 7.2V - 14.4V		
	Voltage	7.2V	
Battery	Charging Time	2Hours or more	
	Running Time 45-50 minutes on full charge		
Packing	Packed in Elega	nt Carrying Case	






Saw Blades

For G3.305 & G3.300UA



Item Code	Length	Width	Thickness	
G3.311A	70mm	24mm	1.0mm	
G3.311B	90mm	24mm	1.0mm	
G3.311C	90mm	13mm	1.0mm	
G3.311D	90mm	19mm	1.0mm	
G3.311E	100mm	22mm	1.0mm	
G3.311F	70mm	14mm	0.5mm	
G3.311G	55mm	15mm	0.5mm	
G3.311H	90mm	12mm	1.2mm	
G3.311J	70mm	25mm	1.0mm	
G3.311K	90mm	25mm	1.2mm	
G3.311L	80mm	20mm	1.2mm	
G3.311M	90mm	22mm	1.2mm	
G3.311N	90mm	19mm	1.2mm	
G3.311P	100mm	22mm	1.2mm	





Saw Blades



Mini Saw Spare Blade For G3.501 & G3.503

Item Code	Length	Width	Thickness	
G3.511A	22mm	08mm	0.5mm	
G3.511B	31mm	9.5mm	0.5mm	
G3.511C	22mm	12mm	0.5mm	
G3.511D	40mm	15mm	0.5mm	





Sternum Spare Blade

For G3.304 (Sternum saw)

Item Code	Length	Width	Thickness
G3.212	38mm	08mm	0.7mm





Reciprocating Spare Blade

For G3.304 (Reciprocating saw)

Item Code	Length	Width	Thickness
G3.312A	78mm	13mm	1.0mm





K-Wire Attachment



K Wire Attachment for G3.302

Item Code	Size (for wire dia)	
G3.219A	0.5mm - 1.5mm	
G3.219F	1.5mm - 3.2mm	







K-Wire Attachment



K Wire Attachment for G3.300UA

Item Code	Size (for wire dia)	
G3.319UG	0.5mm - 1.5mm	
G3.319UF	1.5mm - 3.2mm	









Spares



Chuck with key







Quick Coupling Adapter - for Drill Bits





Reamer Adapter

- for Hexagonal fitting reamer





Reamer Adapter

- for Hudson fitting Reamers







Spares



Battery

For G3.302, G3.305, G3.304, G3.302AR & G3.300UA



Battery Charger

For G3.302, G3.305, G3.304, G3.302AR & G3.300UA



Battery

For G3.500, G3.501 & G3.503,



Battery Charger

For G3.500, G3.501 & G3.503,













Spares



Sterilization Box For Drill/Saw Set

For G3.302, G3.305, G3.304, & G3.302AR







Sterilization Box For Multifunction Machine Set

For G3.300UA







Sterilization Box For Micro Drill/ Micro Saw Set

For G3.500, G3.501







Visit our website for our regularly updated product details:

www.orthocare.in



TERMS & CONDITIONS

- All sizes, dimensions, values etc. in this catalogue are approximate. We reserve the right to upgrade / improve / alter the specifications of our products without notice.
- Use of all our products is restricted to professionals who are qualified and trained to use them. We are not responsible for any damages which may occur due to improper or unauthorized use. In case of malfunctioning of the product, the company's liability is limited to the maximum extent of the cost of equipment.
- Technical specifications are based on test conditions and may vary under actual operating condition

- It is users responsibility to keep the products and packing away from the reach of unauthorised people, children and animals and dispose them off in an environment friendly manner.
- All transactions shall be governed by jurisdiction under Delhi courts.
- You must keep proper record of your sales. This will assist product re-calls (if required).
- All orders are accepted subject to existence of conditions that enable smooth compliance. No liability shall lie against us or any of our representatives, if any orders, remain uncomplied with or partially complied with or delayed compliance for reasons beyond our control.









is simply

S

M

Α

1

S Smooth performance

Meets all needs

A All-round, versatile

R Reasonably priced

Totally reliable









ORIHO CARE

H.O.: 3/10 Kirti Nagar Industrial Area, New Delhi - 110015, INDIA

Phone: \$\infty\$ +919810038246 | +91-11-41424511 | 41424522

E-mail: orthocareindia2009@gmail.com | Web: www.orthocare.in, www.aocare.in

USFDA 510K CLEARED PRODUCTS

(€ CERTIFIED PRODUCTS

ISO 13485 : 2016 CERTIFICATION

LICENCE NO. MFG/MD/2021/000347